

FDA Submission

Claire Vania

Name of the Device: PneumonoTec

Algorithm Description:

1. General Information

Intended Use Statement: Assisting radiologist in detection of pneumonia from Xray

Indications for Use: Intended to use on men and women in the age of 1 to 85 with view position being AP or PA

Device Limitations:

- Prediction time depends upon the processing speed of CPU
- Presence of Infiltration or Atelectasis can affect the prediction

Clinical Impact of Performance: The model had lower precision and higher recall which means that there would be more false positives than false negatives and hence the model should only be used to assist clinicians and not to conclude the presence or absence of a disease

2. Algorithm Design and Function



DICOM Checking Steps:

- Body Part: Chest
- Position: PA or AP
- Modality: DX

Preprocessing Steps:

- Reshaping
- Normalizing
- Repeating the process for all three channels

CNN Architecture: We use a pretrained VGG16 model. The VGG16 model layers were frozen except the last convolution layer and then flattened to pass through fully connected layers

Model: "sequential_1"

Layer (type)	Output Shape	Param #
model_2 (Model)	(None, 7, 7, 512)	14714688
flatten_1 (Flatten)	(None, 25088)	0
dropout_1 (Dropout)	(None, 25088)	0
dense_1 (Dense)	(None, 1024)	25691136
dropout_2 (Dropout)	(None, 1024)	0
dense_2 (Dense)	(None, 512)	524800
dropout_3 (Dropout)	(None, 512)	0
dense_3 (Dense)	(None, 256)	131328
dense_4 (Dense)	(None, 1)	257
Total params: 41,062,209		
Trainable params: 28,707,329		
Non-trainable params: 12,354,880		

3. Algorithm Training

Parameters:

Types of augmentation used during training:

- shear_range = 0.1,
- zoom_range = 0.1,
- rotation_range = 20,
- height_shift_range = 0.1,
- width_shift_range = 0.1

Batch size: 32


Optimizer learning rate: 1e-4

Layers of pre-existing architecture that were frozen: 17 layers of VGG Model


Layers of pre-existing architecture that were fine-tuned: Last two layers

Layers added to pre-existing architecture: flatten, dropout and dense layers as shown in the architecture

Algorithm Training Performance Visualization:



P-R curve




Final Threshold: Threshold was selected where F1 score was maximum

The F1 score from the paper CheXNet: Radiologist-Level Pneumonia Detection on Chest X-Rays with Deep Learning.

	F1 Score (95% CI)
Radiologist 1	0.383 (0.309, 0.453)
Radiologist 2	0.356 (0.282, 0.428)
Radiologist 3	0.365 (0.291, 0.435)
Radiologist 4	0.442 (0.390, 0.492)
Radiologist Avg.	0.387 (0.330, 0.442)
CheXNet	0.435 (0.387, 0.481)

Maximum f1: 0.41352657004830917

Threshold: 0.28675130009651184




The model's F1 score is better than most radiologists and comparable to that of CheXNet.

4. Databases

NIH Chest X-ray Dataset was used which includes 112,000 chest x-rays with disease labels acquired from 30,000 patients.

Description of Training Dataset:

The training set has equal population of pneumonia and non-pneumonia samples with total samples being 2,290. Some training samples.



Description of Validation Dataset:

Validation dataset had 1:4 ratio of pneumonia positive to pneumonia negative cases to represent real world scenario. Total samples are 1430

5. Ground Truth

There are 112,120 X-ray images with disease labels from 30805 unique patients in this dataset. The disease labels were created using Natural Language Processing (NLP) to mine the associated radiological reports. The labels include 14 common thoracic pathologies.

- Consolidation
- Infiltration
- Pneumothorax
- Edema
- Emphysema
- Fibrosis
- Effusion
- Pneumonia
- Pleural Thickening
- Cardiomegaly
- Nodule
- Mass
- Hernia

The biggest limitation of the dataset is that image labels were NLP-extracted so there could be some erroneous labels, but the NLP labeling accuracy is estimated to be >90%

6. FDA Validation Plan

Patient Population Description for FDA Validation Dataset: Patients can be male or female in the age range of 1 – 85 with chest xray being taken in either PA or AP position. Since the dataset should have high correlation between infiltration and Atelectasis with pneumonia, the validation set should not include patients with these diseases.

Ground Truth Acquisition Methodology: The gold standard is by obtaining a lab report for the presence or absence of pneumonia bacteria. Since this is time consuming, we can rely on silver standard which is taking weighted sum of several radiologist defined labels.

Algorithm Performance Standard: The algorithms F1 score should be more than that of an average radiologists(0.387) as mentioned in the ChesXNet paper.